



POSITIONING STRAPS

DISPOSABLE STRAPS FOR POSITIONING
DURING MEDICAL PROCEDURES

PRODUCT CATALOG

PICTURE	NAME	ID	SIZE	PACKAGING	DESCRIPTION	USE
	THIGH STRAPS	1121	10X120CM	200 pcs	SMALL	These thigh straps are used for securing patients during surgery and under anaesthesia. The length of the straps varies to accommodate different sizes of people.
		1140	10X150CM	150 pcs	MEDIUM	
		601100	10X220CM	150 pcs	LARGE	
		1150	10X260CM	100 pcs	X-LARGE	
		1101	10X300CM	100 pcs	XX-LARGE	
	FULL VELCRO LEG STRAP	601102	10X220CM	100 pcs	ONE SIZE	The leg strap with full Velcro is used for securing patients during surgery and under anaesthesia.
	BELLY STRAP	601130	20X220CM	50 pcs	ONE SIZE	This belly strap is used to secure the patient's abdomen during surgery and under anaesthesia.
		1131	25X220CM	50 pcs	ONE SIZE	
	LEG STRAPS	1440	30X160CM	25 pairs	ONE SIZE	These leg straps are used for securing the patient's lower extremity during surgery such as laparoscopic surgery, where the patient is placed on a surgical bed with split leg plates.
	TORSO STRAP	1442	25x200CM	25 pcs	ONE SIZE	Torso straps are used to safely position the arms during ear, nose and throat surgery, when the arms are positioned by the side of the patient.
	ARM/GU STRAP	1112	5X60CM	200 pairs	X-SMALL	Arm/GU straps are used to safely position the patient's arms on the arm table during surgery and under anaesthesia. The straps can also be used for the lower extremity when patients are placed on split surgical beds.
		1111	5X70CM	200 pairs	SMALL	
		601110	5X85CM	200 pairs	MEDIUM	
		1120	5X120CM	200 pairs	LARGE	